

HIGH-ENERGY PET

# SHARPS CONTAINER SHIELD

The Sharps Container Shield for PET is a simple, safe and convenient solution for disposal of syringes that have been contaminated with high-energy radionuclides. The shield is constructed of steel and lined with 1" thick (2.5 cm) lead.

The shield is designed to be used with 039-414 two quart Sharps Containers. It features a lockable sliding cover for container removal and a hinged top door for syringe disposal.

The shield can stand independently or can be recessed into a cabinet or countertop.

## FEATURES:

- *Full 1-inch shielding for PET*
- *Can be built into cabinetry*
- *Lockable cover*



## SPECIFICATIONS:

**Dimensions:** 12" h x 8.75" dia (30.5 x 22 cm)

**Lead Shielding:** 1" thick (2.5 cm)

**Security:** Key-locked

**Weight:** 175 lb (79.4 kg)

**Shipping Weight:** 222 lb (100.6 kg)

**039-412** Sharps Container Shield, PET, 1" lead .....\$1,450.00  
*Uses one 039-413 Sharps Container*

## Related:

**039-413** Sharps Container, 2 qt., 30/pkg ..... 175.00  
*Fits 039-412*

*Prices subject to change without notice.*

[www.jgravengaard.com](http://www.jgravengaard.com) 2007-10-08

**J. Gravengaard Corporation**

Medical Imaging Sales & Service

PO Box 20171 Portland, Oregon 97294 Nationwide (800)796-SCAN

Local (503)255-0535 Fax (503)255-7047