

INJECTION / RESTING CHAIR

A comfortable rest between injection and imaging



Most protocols require that prior to imaging, a patient rest for up to one hour after an injection of FDG F-18. Typical injection chairs are not designed for comfort, making rest difficult.

The Injection/Resting Chair is a comfortable resting spot for patients to wait between injection and imaging. The chair achieves three recline positions and features an infinitely adjustable back. For patient safety, deep recline and Trendelenburg positions can only be achieved by the attendant. A quick release handle that spans the entire width of chair back allows effortless position change from either side of the chair.

SPECIFICATIONS:

Dimensions: 35" w x 49" h (89 x 125 cm)

Seat: 20" depth x 25" w x 21" h (51x64x53cm)

Back: 25" w x 34" h, above seat (64 x 86 cm)

Upholstery: Blueridge, fire retardant

Frame: Welded powder coated steel

Casters: Four 5" heavy-duty with positive locking swivel and wheel brakes

Warranty: 5-year on frame

Weight Capacity: 450 lb (204 kg)

Weight: 111 lb (50.4 kg)

214-210 Chair, Injection / Resting \$1,495.00

Prices subject to change without notice.

www.jgravengaard.com 2007-07-09

J. Gravengaard Corporation

Medical Imaging Sales & Service

PO Box 20171 Portland, Oregon 97294 Nationwide (800)796-SCAN Local (503)255-0535 Fax (503)255-7047